



Plastic Surgery Instruments Catalog

www.healermedinstruments.com
Email: Healermedins@gmail.com
Whatsapp & Call: +92 3034515151



Let Us Know if you could not find what you are looking for

CONTENTS

1 Scalpel Handles	16 Rhinoplasty Raspatories
2 Breast Areola Markers	17 Nasal Retractors
3 Measuring Caliper	18 Abdominal Retractors
4 Dematomes	19 Fiber Optic Retractors
5 Plastic Surgery Scissors	20 Breast Dissectors
6 Hemastatic Forceps	21 Forehead Lift Dissectors
7 Facelift Forceps	22 Bone Mallets
8 Needle Forceps	23 Cartilage Instruments
9 Micro Scissor/ Forceps	
10 Tissue Forceps	
11 Skin Hooks	
12 Chisels & Osteotomes	
13 Septum & Periosteal Elevators	
14 Rhinoplasty Knives	
15 Nasal Speculums/ Saws	



HMI-PS-101

No. 3



HMI-PS-102

No. 4



HMI-PS-103

No. 3R



HMI-PS-104

No. 3RG



HMI-PS-105

No. 7K



HMI-PS-106

No. 7



Solz Aerola Marker

HMI-PS-201	32mm
HMI-PS-202	38mm
HMI-PS-203	40mm
HMI-PS-204	42mm



Freeman Aerola Marker Double
End

HMI-PS-205	38/42mm
HMI-PS-206	42/45mm



Freeman Aerola
Marker

HMI-PS-207
HMI-PS-208
HMI-PS-209
HMI-PS-210
HMI-PS-211
HMI-PS-212
HMI-PS-213
HMI-PS-214



Grossman Aerola marker
Four sizes in One

HMI-PS-215	40/44/48/50mm
------------	---------------



McKissock Keyhole

HMI-PS-216	38mm
HMI-PS-217	42mm
HMI-PS-218	45mm



HMI-PS-301

Precision Caliper 10cm



HMI-PS-302

Jameson Caliper 10cm




HMI-PS-304

Thrope Caliper 11cm, 80mm



HMI-PS-303

Castroviejo Caliper, 9cm,
40mm



HMI-PS-305

S/ steel Ruler 12"

HMI-PS-306

S/ steel Ruler 6"



HMI-PS-401
Tiersch 20cm



HMI-PS-402
Schink 30cm



HMI-PS-403
Humby 31cm



HMI-PS-404
Watson 31cm



HMI-PS-405
Silver 19cm



HMI-PS-406
Schink Plate



HMI-PS-407
Disposable Blade

SKIN GRAFT HUMBY HANDLE SPECIFICATIONS

Stainless steel handle Autoclavable

Length of 28-33cm

Compatible with blades of 15-16cm

Fixed polished (non movable) steel guard

Adjustable steel rod with lock at opposite end

Knurled control knob on the adjustable rod that sets depth

Skin thickness range 0.00mm-1mm with increment intervals of 0.2-0.3mm

Locking leaf blade holder that opens at $>60^\circ$

Safety catch on the leaf

1-2 year warranty

Packed in Wooden box



HMI-PS-408



Standard

Size	Straight	Curved
11.5 cm	HMI-PS-501	HMI-PS-507
14 cm	HMI-PS-502	HMI-PS-508
15 cm	HMI-PS-503	HMI-PS-509
16 cm	HMI-PS-504	HMI-PS-510
18 cm	HMI-PS-505	HMI-PS-511
20 cm	HMI-PS-506	HMI-PS-512

TC

Size	Straight	Curved
11.5cm	HMI-PS-513	HMI-PS-519
14 cm	HMI-PS-514	HMI-PS-520
15cm	HMI-PS-515	HMI-PS-521
16cm	HMI-PS-516	HMI-PS-522
18 cm	HMI-PS-517	HMI-PS-523
20cm	HMI-PS-518	HMI-PS-524

SUPER CUT TC

Size	Straight	Curved
11.5 cm	HMI-PS-525	HMI-PS-531
14 cm	HMI-PS-526	HMI-PS-532
15 cm	HMI-PS-527	HMI-PS-533
16 cm	HMI-PS-528	HMI-PS-534
18 cm	HMI-PS-529	HMI-PS-535
20 cm	HMI-PS-530	HMI-PS-536

Supercut

Size	Straight	Curved
11.5 cm	HMI-PS-537	HMI-PS-542
14 cm	HMI-PS-538	HMI-PS-543
16 cm	HMI-PS-539	HMI-PS-544
18 cm	HMI-PS-540	HMI-PS-545
20 cm	HMI-PS-541	HMI-PS-546

Let Us know if you could not find what you are looking for



Standard

Size	Straight	Curved
11.5cm	HMI-PS-547	HMI-PS-553
14 cm	HMI-PS-548	HMI-PS-554
15cm	HMI-PS-549	HMI-PS-555
16cm	HMI-PS-550	HMI-PS-556
18 cm	HMI-PS-551	HMI-PS-557
20cm	HMI-PS-552	HMI-PS-558

TC

Size	Straight	Curved
11.5 cm	HMI-PS-560	HMI-PS-566
14 cm	HMI-PS-561	HMI-PS-567
15 cm	HMI-PS-562	HMI-PS-568
16 cm	HMI-PS-563	HMI-PS-569
18 cm	HMI-PS-564	HMI-PS-570
20 cm	HMI-PS-565	HMI-PS-571

SUPER
CUT

TC

Size	Straight	Curved
11.5 cm	HMI-PS-572	HMI-PS-577
14 cm	HMI-PS-573	HMI-PS-578
16 cm	HMI-PS-574	HMI-PS-579
18 cm	HMI-PS-575	HMI-PS-580
20 cm	HMI-PS-576	HMI-PS-581

Supercut

Size	Straight	Curved
11.5 cm	HMI-PS-582	HMI-PS-587
14 cm	HMI-PS-583a	HMI-PS-588a
15 cm	HMI-PS-583	HMI-PS-588
16 cm	HMI-PS-584	HMI-PS-589
18 cm	HMI-PS-585	HMI-PS-590
20 cm	HMI-PS-586	HMI-PS-591



Let Us know if you could not find what you are looking for

Standard



Size	Straight	Curved
11.5 cm	HMI-PS-592	HMI-PS-598
14 cm	HMI-PS-593	HMI-PS-599
15 cm	HMI-PS-594	HMI-PS-600
16 cm	HMI-PS-595	HMI-PS-601
18 cm	HMI-PS-596	HMI-PS-602
20 cm	HMI-PS-597	HMI-PS-603

TC



Size	Straight	Curved
11.5cm	HMI-PS-604	HMI-PS-610
14 cm	HMI-PS-605	HMI-PS-611
15cm	HMI-PS-606	HMI-PS-612
16cm	HMI-PS-607	HMI-PS-613
18 cm	HMI-PS-608	HMI-PS-614
20cm	HMI-PS-609	HMI-PS-615

SUPER CUT

TC



Size	Straight	Curved
11.5 cm	HMI-PS-616	HMI-PS-622
14 cm	HMI-PS-617	HMI-PS-623
15 cm	HMI-PS-618	HMI-PS-624
16 cm	HMI-PS-619	HMI-PS-625
18 cm	HMI-PS-620	HMI-PS-626
20 cm	HMI-PS-621	HMI-PS-627

Supercut



Size	Straight	Curved
11.5 cm	HMI-PS-628	HMI-PS-633
14 cm	HMI-PS-629	HMI-PS-634
15 cm	HMI-PS-629a	HMI-PS-634a
16 cm	HMI-PS-630	HMI-PS-635
18 cm	HMI-PS-631	HMI-PS-636
20 cm	HMI-PS-632	HMI-PS-637

Let Us know if you could not find what you are looking for

Standard



Size	Straight	Curved
14 cm	HMI-PS-001	HMI-PS-006
15 cm	HMI-PS-002	HMI-PS-007
16 cm	HMI-PS-003	HMI-PS-008
18 cm	HMI-PS-004	HMI-PS-009
20 cm	HMI-PS-005	HMI-PS-010

TC



Size	Straight	Curved
14 cm	HMI-PS-011	HMI-PS-016
15 cm	HMI-PS-012	HMI-PS-017
16 cm	HMI-PS-013	HMI-PS-018
18 cm	HMI-PS-014	HMI-PS-019
20 cm	HMI-PS-015	HMI-PS-020

SUPER
CUT

TC



Size	Straight	Curved
14 cm	HMI-PS-021	HMI-PS-026
15 cm	HMI-PS-022	HMI-PS-027
16 cm	HMI-PS-023	HMI-PS-028
18 cm	HMI-PS-024	HMI-PS-029
20 cm	HMI-PS-025	HMI-PS-030

Supercut



Size	Straight	Curved
14 cm	HMI-PS-031	HMI-PS-036
15 cm	HMI-PS-032	HMI-PS-037
16 cm	HMI-PS-033	HMI-PS-038
18 cm	HMI-PS-034	HMI-PS-039
20 cm	HMI-PS-035	HMI-PS-040

Let Us know if you could not find what you are looking for

Standard

Size	Straight	Curved
15 cm	HMI-PS-041	HMI-PS-044
16 cm	HMI-PS-042	HMI-PS-045
18 cm	HMI-PS-043	HMI-PS-046

TC

Size	Straight	Curved
15 cm	HMI-PS-047	HMI-PS-050
16 cm	HMI-PS-048	HMI-PS-051
18 cm	HMI-PS-049	HMI-PS-052

SUPER
CUT TC

Size	Straight	Curved
15 cm	HMI-PS-053	HMI-PS-056
16 cm	HMI-PS-054	HMI-PS-057
18 cm	HMI-PS-055	HMI-PS-058

Supercut

Size	Straight	Curved
15 cm	HMI-PS-059	HMI-PS-062
16 cm	HMI-PS-060	HMI-PS-063
18 cm	HMI-PS-061	HMI-PS-064



Let Us know if you could not find what you are looking for



Standard

Size	Straight	Curved
15 cm	HMI-PS-638	HMI-PS-640
16 cm	HMI-PS-639	HMI-PS-641



TC

Size	Straight	Curved
15 cm	HMI-PS-642	HMI-PS-644
16 cm	HMI-PS-643	HMI-PS-645



SUPER CUT TC

Size	Straight	Curved
15 cm	HMI-PS-646	HMI-PS-648
16 cm	HMI-PS-647	HMI-PS-649



Supercut

Size	Straight	Curved
15 cm	HMI-PS-650	HMI-PS-652
16 cm	HMI-PS-651	HMI-PS-653

Let Us know if you could not find what you are looking for



Standard

Size	Straight	Curved
15 cm	HMI-PS-654	HMI-PS-656
16 cm	HMI-PS-655	HMI-PS-657



TC

Size	Straight	Curved
10.5 cm	HMI-PS-658	HMI-PS-660
11.5 cm	HMI-PS-659	HMI-PS-661



Supercut

Size	Straight	Curved
10.5 cm	HMI-PS-662	HMI-PS-664
11.5 cm	HMI-PS-663	HMI-PS-665

SUPER CUT TC

Size	Straight	Curved
10.5 cm	HMI-PS-662	HMI-PS-664
11.5 cm	HMI-PS-663	HMI-PS-665

Let Us know if you could not find what you are looking for



Standard

Size	Straight	Curved
14 cm	HMI-PS-670	HMI-PS-675
15 cm	HMI-PS-671	HMI-PS-676
18 cm	HMI-PS-672	HMI-PS-677
20 cm	HMI-PS-673	HMI-PS-678
23 cm	HMI-PS-674	HMI-PS-679

TC

Size	Straight	Curved
14 cm	HMI-PS-680	HMI-PS-685
15 cm	HMI-PS-681	HMI-PS-686
18 cm	HMI-PS-682	HMI-PS-687
20 cm	HMI-PS-683	HMI-PS-688
23 cm	HMI-PS-684	HMI-PS-689

SUPER CUT TC

Size	Straight	Curved
14 cm	HMI-PS-690	HMI-PS-695
15 cm	HMI-PS-691	HMI-PS-696
18 cm	HMI-PS-692	HMI-PS-697
20 cm	HMI-PS-693	HMI-PS-698
23 cm	HMI-PS-694	HMI-PS-699

Supercut

Size	Straight	Curved
14 cm	HMI-PS-700	HMI-PS-705
15 cm	HMI-PS-701	HMI-PS-706
18 cm	HMI-PS-702	HMI-PS-707
20 cm	HMI-PS-703	HMI-PS-708
23 cm	HMI-PS-704	HMI-PS-709

Let Us know if you could not find what you are looking for



Standard

Size	Straight	Curved
14 cm	HMI-PS-710	HMI-PS-715
15 cm	HMI-PS-711	HMI-PS-716
18 cm	HMI-PS-712	HMI-PS-717
20 cm	HMI-PS-713	HMI-PS-718
23 cm	HMI-PS-714	HMI-PS-719

TC

Size	Straight	Curved
14 cm	HMI-PS-720	HMI-PS-725
15 cm	HMI-PS-721	HMI-PS-726
18 cm	HMI-PS-722	HMI-PS-727
20 cm	HMI-PS-723	HMI-PS-728
23 cm	HMI-PS-724	HMI-PS-729

SUPER CUT TC

Size	Straight	Curved
14 cm	HMI-PS-730	HMI-PS-735
15 cm	HMI-PS-731	HMI-PS-736
18 cm	HMI-PS-732	HMI-PS-737
20 cm	HMI-PS-733	HMI-PS-738
23 cm	HMI-PS-734	HMI-PS-739

Supercut

Size	Straight	Curved
14 cm	HMI-PS-740	HMI-PS-745
15 cm	HMI-PS-741	HMI-PS-746
18 cm	HMI-PS-742	HMI-PS-747
20 cm	HMI-PS-743	HMI-PS-748
23 cm	HMI-PS-744	HMI-PS-749

Let Us know if you could not find what you are looking for



Standard

Size	Straight	Curved
11.5 cm	HMI-PS-750	HMI-PS-756
14 cm	HMI-PS-751	HMI-PS-757
15 cm	HMI-PS-752	HMI-PS-758
16 cm	HMI-PS-753	HMI-PS-759
18 cm	HMI-PS-754	HMI-PS-760
20 cm	HMI-PS-755	HMI-PS-761

TC

Size	Straight	Curved
11.5cm	HMI-PS-762	HMI-PS-768
14 cm	HMI-PS-763	HMI-PS-769
15cm	HMI-PS-764	HMI-PS-770
16cm	HMI-PS-765	HMI-PS-771
18 cm	HMI-PS-766	HMI-PS-772
20cm	HMI-PS-767	HMI-PS-773

SUPER CUT TC

Size	Straight	Curved
11.5 cm	HMI-PS-774	HMI-PS-780
14 cm	HMI-PS-775	HMI-PS-781
15 cm	HMI-PS-776	HMI-PS-782
16 cm	HMI-PS-777	HMI-PS-783
18 cm	HMI-PS-778	HMI-PS-784
20 cm	HMI-PS-779	HMI-PS-785

Supercut

Size	Straight	Curved
11.5 cm	HMI-PS-786	HMI-PS-791
14 cm	HMI-PS-787	HMI-PS-792
16 cm	HMI-PS-788	HMI-PS-793
18 cm	HMI-PS-789	HMI-PS-794
20 cm	HMI-PS-790	HMI-PS-795

Let Us know if you could not find what you are looking for



Standard

Size	Straight	Curved
11.5 cm	HMI-PS-796	HMI-PS-797



TC

Size	Straight	Curved
11.5 cm	HMI-PS-798	HMI-PS-799

SUPER CUT TC


Size	Straight	Curved
11.5 cm	HMI-PS-780	HMI-PS-781



Supercut

Size	Straight	Curved
11.5 cm	HMI-PS-782	HMI-PS-783

Let Us know if you could not find what you are looking for



Standard

Size	Straight	Curved
15.5 cm	HMI-PS-784	HMI-PS-786
18 cm	HMI-PS-785	HMI-PS-787



TC

Size	Straight	Curved
15.5 cm	HMI-PS-788	HMI-PS-790
18 cm	HMI-PS-789	HMI-PS-791




SuperCut

Size	Straight	Curved
15.5 cm	HMI-PS-792	HMI-PS-794
18 cm	HMI-PS-793	HMI-PS-795

SUPER CUT TC

Size	Straight	Curved
15.5 cm	HMI-PS-792	HMI-PS-794
18 cm	HMI-PS-793	HMI-PS-795

Let Us know if you could not find what you are looking for



Standard

Size	Straight	Curved
14.5 cm	HMI-PS-780	HMI-PS-781

TC

Size	Straight	Curved
14.5 cm	HMI-PS-782	HMI-PS-783



SUPER CUT TC


Size	Straight	Curved
14.5 cm	HMI-PS-784	HMI-PS-785

Supercut

Size	Straight	Curved
14.5 cm	HMI-PS-786	HMI-PS-787



Let Us know if you could not find what you are looking for



Standard

Size	Straight	Curved
14 cm	HMI-PS-788	HMI-PS-789



TC

Size	Straight	Curved
14 cm	HMI-PS-790	HMI-PS-791

SUPER CUT TC

Size	Straight	Curved
14 cm	HMI-PS-792	HMI-PS-793

Supercut

Size	Straight	Curved
14 cm	HMI-PS-794	HMI-PS-795

Let Us know if you could not find what you are looking for



Standard

Size	Straight	Curved
16 cm	HMI-PS-796	HMI-PS-797



TC

Size	Straight	Curved
16 cm	HMI-PS-798	HMI-PS-799


SUPER CUT TC

Size	Straight	Curved
16 cm	HMI-PS-800	HMI-PS-801

Supercut

Size	Straight	Curved
16 cm	HMI-PS-802	HMI-PS-803

Let Us know if you could not find what you are looking for



Standard

Size	Straight	Curved
11.5 cm	HMI-PS-804	HMI-PS-806
13 cm	HMI-PS-805	HMI-PS-807

TC

Size	Straight	Curved
11.5 cm	HMI-PS-808	HMI-PS-810
13 cm	HMI-PS-809	HMI-PS-811



SUPER CUT TC

Size	Straight	Curved
11.5 cm	HMI-PS-812	HMI-PS-814
13 cm	HMI-PS-813	HMI-PS-815

Supercut



Size	Straight	Curved
11.5 cm	HMI-PS-816	HMI-PS-818
13 cm	HMI-PS-817	HMI-PS-819

Let Us know if you could not find what you are looking for

Standard

Size	Straight	Curved
9.5 cm	HMI-PS-820	HMI-PS-821



TC

Size	Straight	Curved
9.5 cm	HMI-PS-822	HMI-PS-823



SUPER CUT TC

Size	Straight	Curved
9.5 cm	HMI-PS-824	HMI-PS-825

Supercut

Size	Straight	Curved
9.5 cm	HMI-PS-826	HMI-PS-827



Let Us know if you could not find what you are looking for



Standard

Size	Straight	Curved
14.5 cm	HMI-PS-828	HMI-PS-829



TC

Size	Straight	Curved
14.5 cm	HMI-PS-830	HMI-PS-831



SUPER CUT TC

Size	Straight	Curved
14.5 cm	HMI-PS-832	HMI-PS-833

Supercut

Size	Straight	Curved
14.5 cm	HMI-PS-834	HMI-PS-835

Let Us know if you could not find what you are looking for



Standard

Size	Straight	Curved
11.5 cm	HMI-PS-836	HMI-PS-837



TC

Size	Straight	Curved
11.5 cm	HMI-PS-838	HMI-PS-839

SUPER CUT TC


Size	Straight	Curved
11.5 cm	HMI-PS-840	HMI-PS-841



Supercut

Size	Straight	Curved
11.5 cm	HMI-PS-842	HMI-PS-843

Let Us know if you could not find what you are looking for



Standard

Size	Straight	Curved
12 cm	HMI-PS-844	HMI-PS-846
15 cm	HMI-PS-845	HMI-PS-847

TC

Size	Straight	Curved
12 cm	HMI-PS-848	HMI-PS-850
15 cm	HMI-PS-849	HMI-PS-851



SUPER
CUT TC

Size	Straight	Curved
12 cm	HMI-PS-852	HMI-PS-854
15 cm	HMI-PS-853	HMI-PS-855

Supercut

Size	Straight	Curved
12 cm	HMI-PS-856	HMI-PS-858
15 cm	HMI-PS-857	HMI-PS-859



Standard

Size	Straight	Curved
12 cm	HMI-PS-860	HMI-PS-862
18 cm	HMI-PS-861	HMI-PS-863



TC

Size	Straight	Curved
12 cm	HMI-PS-864	HMI-PS-866
18 cm	HMI-PS-865	HMI-PS-867

SUPER CUT TC




Supercut

Size	Straight	Curved
12 cm	HMI-PS-868	HMI-PS-870
18 cm	HMI-PS-869	HMI-PS-871

Size	Straight	Curved
12 cm	HMI-PS-872	HMI-PS-874
18 cm	HMI-PS-873	HMI-PS-875

Let Us know if you could not find what you are looking for



Size	Straight	Curved
18 cm	HMI-PS-876	HMI-PS-877



Size	Straight	Curved
18 cm	HMI-PS-878	HMI-PS-879

Size	Straight	Curved
18 cm	HMI-PS-880	HMI-PS-881

Size	Straight	Curved
18 cm	HMI-PS-882	HMI-PS-883

Let Us know if you could not find what you are looking for



Standard

Size	Straight	Curved
14 cm	HMI-PS-884	HMI-PS-885



TC

Size	Straight	Curved
14 cm	HMI-PS-886	HMI-PS-887



SUPER CUT TC

Size	Straight	Curved
14 cm	HMI-PS-888	HMI-PS-889



Supercut

Size	Straight	Curved
14 cm	HMI-PS-890	HMI-PS-891

Let Us know if you could not find what you are looking for



Standard

Size	Straight	Curved
14.5 cm	HMI-PS-890	HMI-PS-892
18 cm	HMI-PS-891	HMI-PS-893



TC

Size	Straight	Curved
14.5 cm	HMI-PS-894	HMI-PS-896
18 cm	HMI-PS-895	HMI-PS-897



SUPER CUT TC

Size	Straight	Curved
14.5 cm	HMI-PS-898	HMI-PS-900
18 cm	HMI-PS-899	HMI-PS-901



Supercut

Size	Straight	Curved
14.5 cm	HMI-PS-902	HMI-PS-904
18 cm	HMI-PS-903	HMI-PS-905

Let Us know if you could not find what you are looking for



Standard

Size	Straight	Curved
13 cm	HMI-PS-906	HMI-PS-907



TC

Size	Straight	Curved
13 cm	HMI-PS-908	HMI-PS-909



SUPER CUT TC

Size	Straight	Curved
13 cm	HMI-PS-910	HMI-PS-911

Supercut

Size	Straight	Curved
13 cm	HMI-PS-912	HMI-PS-913

Let Us know if you could not find what you are looking for



Standard

Size	Straight	Curved
13 cm	HMI-PS-914	HMI-PS-915

TC

Size	Straight	Curved
13 cm	HMI-PS-916	HMI-PS-917



SUPER CUT TC

Size	Straight	Curved
13 cm	HMI-PS-918	HMI-PS-919



Supercut

Size	Straight	Curved
13 cm	HMI-PS-920	HMI-PS-921

Let Us know if you could not find what you are looking for



Standard

Size	Straight	Curved
18 cm	HMI-PS-920	HMI-PS-921



TC

Size	Straight	Curved
18 cm	HMI-PS-922	HMI-PS-923

SUPER CUT TC

Size	Straight	Curved
18 cm	HMI-PS-924	HMI-PS-925

Supercut

Size	Straight	Curved
18 cm	HMI-PS-926	HMI-PS-927



Let Us know if you could not find what you are looking for



Standard

Size	Straight	Curved
13.5 cm	HMI-PS-928	HMI-PS-929



TC

Size	Straight	Curved
13.5 cm	HMI-PS-930	HMI-PS-931

SUPER CUT TC

Size	Straight	Curved
13.5 cm	HMI-PS-932	HMI-PS-933

Supercut

Size	Straight	Curved
13.5 cm	HMI-PS-934	HMI-PS-935



Let Us know if you could not find what you are looking for



Standard

Size	Straight	Curved
19 cm	HMI-PS-936	HMI-PS-937



TC

Size	Straight	Curved
19 cm	HMI-PS-938	HMI-PS-939



SUPER CUT TC

Size	Straight	Curved
19 cm	HMI-PS-940	HMI-PS-941

Supercut

Size	Straight	Curved
19 cm	HMI-PS-942	HMI-PS-943

Let Us know if you could not find what you are looking for



HARTMANN

MICRO
MOSQUITO



Size	Straight	Curved
10 cm	HMI-PS-1001	HMI-PS-1002

Size	Straight	Curved
10 cm	HMI-PS-1003	HMI-PS-1005
12 cm	HMI-PS-1004	HMI-PS-1006



Halsted Mosquito

Crile



Size	Straight	Curved
12.5 cm	HMI-PS-1007	HMI-PS-1008

Size	Straight	Curved
14 cm	HMI-PS-1009	HMI-PS-1011
16 cm	HMI-PS-1010	HMI-PS-1012



Allis tiusse

Size	Straight
15 cm	HMI-PS-1013
18 cm	HMI-PS-1014

Babcock



Size	Straight
16 cm	HMI-PS-1015
18 cm	HMI-PS-1016
20 cm	HMI-PS-1017



BACKHAUS

Size	Straight
9 cm	HMI-PS-1018
11 cm	HMI-PS-1019
13 cm	HMI-PS-1020

Baby Mixter




Size	Straight
14 cm	HMI-PS-1021
18 cm	HMI-PS-1022



D'Assumpcao Face Lift Forceps

Size	Straight
16 cm	HMI-PS-1101



Pitanguy Forceps

Size	Straight
15 cm	HMI-PS-1102
26 cm	HMI-PS-1103



Marten Flap Forceps

Size	Straight
15 cm	HMI-PS-1104

Pitanguy Flap Grasping Forceps



Size	Straight
15 cm	HMI-PS-1105

Let Us know if you could not find what you are looking for

Needle Holders

DERF



Size	TC	Standard
12 cm	HMI-PS-1201	HMI-PS-1202

HALSEY



Size	TC	Standard
13 cm	HMI-PS-1203	HMI-PS-1204

BAUMGARTNER



Size	TC	Standard
14 cm	HMI-PS-1205	HMI-PS-1206

WEBSTER



Size	TC	Standard
13 cm	HMI-PS-1207	HMI-PS-1208

038A

Crile Wood



Size	TC	Standard
15 cm	HMI-PS-1209	HMI-PS-1213
18 cm	HMI-PS-1210	HMI-PS-1214
20 cm	HMI-PS-1211	HMI-PS-1215
22 cm	HMI-PS-1212	HMI-PS-1216

Mayo Hegar



Size	TC	Standard
14 cm	HMI-PS-1217	HMI-PS-1222
16 cm	HMI-PS-1218	HMI-PS-1223
18 cm	HMI-PS-1219	HMI-PS-1224
20 cm	HMI-PS-1220	HMI-PS-1225
24 cm	HMI-PS-1221	HMI-PS-1225

Gillies



Size	TC	Standard
16 cm	HMI-PS-1225	HMI-PS-1226

Olsen Hegar



Size	TC	Standard
14 cm	HMI-PS-1227	HMI-PS-1229
18 cm	HMI-PS-1228	HMI-PS-1230

Let Us know if you could not find what you are looking for



Micro Vascular

TC

Size	Smooth	Serrated
14 cm	HMI-PS-1231	HMI-PS-1233
18 cm	HMI-PS-1232	HMI-PS-1233
20 cm	HMI-PS-1233	HMI-PS-1234

Ryder



TC

Size	TC
13 cm	HMI-PS-1235
15 cm	HMI-PS-1236
18 cm	HMI-PS-1237
20 cm	HMI-PS-1238
24 cm	HMI-PS-1239



Sarat

TC



Standard

Size	TC	Standard
18 cm	HMI-PS-1240	HMI-PS-1241

Nievert



TC

Size	Straight
13 cm	HMI-PS-1242



Size	Smooth	Serrated
14 cm	HMI-PS-1243	HMI-PS-1246
18 cm	HMI-PS-1244	HMI-PS-1247
20 cm	HMI-PS-1245	HMI-PS-1248



Size	TC	Standard
16 cm	HMI-PS-1249	HMI-PS-1252
18 cm	HMI-PS-1250	HMI-PS-1253
20 cm	HMI-PS-1251	HMI-PS-1254



Size	Straight
18 cm	HMI-PS-1255
20 cm	HMI-PS-1256
24 cm	HMI-PS-1257



Size	Straight	Curved
14 cm	HMI-PS-1258	HMI-PS-1259



Micro Scissor 15cm

Standard



Micro Blades
Str, 15cm
HMI-PS-1301



Micro Blades
Cvd, 15cm
HMI-PS-1302

TC



Micro Blades
Str, 15cm
HMI-PS-1303



Micro Blades
Cvd, 15cm
HMI-PS-1304



Micro Needle Holder 15cm

Standard

Micro Tip
Str, 15cm
HMI-PS-1305

Micro Blades
Cvd, 15cm
HMI-PS-1306

TC

Micro Blades
Str, 15cm
HMI-PS-1307

Micro Blades
Cvd, 15cm
HMI-PS-1308

Let Us know if you could not find what you are looking for



Adson

Adson
With teeth



Size	Standard	TC
12 cm	HMI-PS-1401	HMI-PS-1403
15 cm	HMI-PS-1402	HMI-PS-1404

Size	Standard	TC
10 cm	HMI-PS-1405	HMI-PS-1407
15 cm	HMI-PS-1406	HMI-PS-1408



Micro Adson



Adson Brown
7x7



Size	W/O Teeth	With teeth
12 cm	HMI-PS-1409	HMI-PS-1410

Size	Standard	TC
12 cm	HMI-PS-141	HMI-PS-1412



Gillies

TC Gillies



Size	W/O Teeth	W/ teeth
15 cm	HMI-PS-1413	HMI-PS-1414

Size	W/O Teeth	W/ Teeth
15 cm	HMI-PS-1415	HMI-PS-1416



McIndoe

TC McIndoe



Size	W/O Teeth	With teeth
15 cm	HMI-PS-1417	HMI-PS-1418

Size	W/O Teeth	W/ Teeth
15 cm	HMI-PS-1419	HMI-PS-1420

**Potts
SMith**



Size	TC
15 cm	HMI-PS-1421
18 cm	HMI-PS-1422

Rees Bayonet



Size	TC
16 cm	HMI-PS-1423

**Brown
adson
angular**



Size	TC
12 cm	HMI-PS-1424

**Grazer Blepharoplasty
Forceps**



Size	TC
12 cm	HMI-PS-1425

Debakey



Size	1.5mm tip	2.0mm tip
16 cm	HMI-PS-1426	HMI-PS-1426

**Castrovijo Suture
Forceps**



W/O Teeth	W Teeth
HMI-PS-1427	HMI-PS-1428

Let Us know if you could not find what you are looking for



GILLIES SKIN HOOK
18cm
3mm

HMI-PS-1501

KLIENERT KUT
Skin Hook 16cm
Sharp

HMI-PS-1502

Joseph Hook
Sharp Prongs

HMI-PS-1503	2mm
HMI-PS-1504	5mm
HMI-PS-1505	7mm
HMI-PS-1506	10mm

Rollet Retractor 4 Prongs
Sharp 14cm

HMI-PS-1507

Senn-Miller
16cm

HMI-PS-1508

HMI-PS-1509

Sharp

Blunt



Freeman rake Retractor
10cm sharp prongs

HMI-PS-1510



Kilner Double Nasal
Hook

HMI-PS-1511



MARTINEZ KILNR
Retractor 11cm
10mm wide

HMI-PS-1512

Let Us know if you could not find what you are looking for




Adson hook 20cm
Sharp

HMI-PS-1513



Bergman Hook Sharp
14cm

HMI-PS-1514




Gillies Dura
Skin Hook

HMI-PS-1515	3mm
HMI-PS-1516	8mm



MEYERDING RETRAC-
TOR18CM

HMI-PS-1517



Humby Swivel
Hook Sharp

HMI-PS-1518



Freeman rake Retractor
10cm sharp prongs

HMI-PS-1519



MARTINEZ KILNR
Retractor 11cm
10mm wide

HMI-PS-1520

Let Us know if you could not find what you are looking for



ALM Retractor 7 cm 4x4
teeth

HMI-PS-1521



Yancoskie 18cm
150mm

HMI-PS-1522



Freeman 20cm

HMI-PS-1523	12mm
HMI-PS-1524	25 mm



Ferreira 4"
x=23mm
y=120mm

HMI-PS-1525



Weitlaner Retractor

HMI-PS-1526	10cm Sharp
HMI-PS-1527	13cm Sharp
HMI-PS-1528	10cm Blunt
HMI-PS-1529	13cm Blunt

Let Us know if you could not find what you are looking for



HMI-PS-1601



HMI-PS-1602



HMI-PS-1603

4mm

7mm

12mm

HMI-PS-1604



HMI-PS-1605

Killian Claus 16cm
5mm

Killian 16cm
5mm

Cottle 15cm

Cottle 16cm, 6mm

HMI-PS-1606



Cottle 16cm, 4mm

HMI-PS-1607





Cottle 18cm

4mm



HMI-PS-065

7mm



HMI-PS-066

9mm



HMI-PS-067

12mm



HMI-PS-068



Cottle 18cm

4mm



HMI-PS-070

6mm



HMI-PS-071

6mm



HMI-PS-072

7mm



HMI-PS-073

8mm




HMI-PS-074

Healermmed instruments ISO 13480 - CE

Let Us know if you could not find what you are looking for




10mm




HMI-PS-1608

12mm




HMI-PS-1609

14mm



HMI-PS-1610

16mm



HMI-PS-1611

Cinelli 18cm



10mm



HMI-PS-1612

12mm




HMI-PS-1613

14mm




HMI-PS-1614

Rubin 18cm




HMI-PS-1615
Straight


Silver 18cm




HMI-PS-1616
Right



HMI-PS-1617
Left



HMI-PS-1618
Right



HMI-PS-1619
Left

Neivert Anderson 20cm



Neivert 20cm
HMI-PS-1620



Cottle 18cm
HMI-PS-1701



**Joseph 16cm
4mm**
HMI-PS-1702



Hale 16cm
HMI-PS-1703




**Ballenger-Hajek
18cm**
HMI-PS-1704



**Freer Blunt
18cm**
HMI-PS-1705



**Freer Blunt
18cm**
HMI-PS-1706



Cottle 18cm

HMI-PS-1707



Cottle 20cm

HMI-PS-1708




Masing 20cm

HMI-PS-1709






Obwegeser
17.5cm



6mm
HMI-PS-1801




8mm
HMI-PS-1802



10mm
HMI-PS-1803



12mm
HMI-PS-1804




Obwegeser
17.5cm
HMI-PS-1805



Obwegeser
18cm
HMI-PS-1806



Obwegeser
22cm
HMI-PS-1807



Obwegeser
17.5cm
HMI-PS-1808



Obwegeser
17.5cm
HMI-PS-1809
7mm
HMI-PS-1810
10mm



Obwegeser
20cm
HMI-PS-1811



5mm
HMI-PS-1812

Obwegeser
20cm



8mm
HMI-PS-1813



10mm
HMI-PS-1814



Obwegeser
18cm
HMI-PS-1816



HMI-PS-1815



HMI-PS-1901
Converse button
end Knife 16cm Curved



HMI-PS-1902
Cottle Septum
Knife 14cm



HMI-PS-1903
Cottle Septum
Knife 14cm



HMI-PS-1904
Joseph Angled
button end Knife 17.5cm




HMI-PS-1905
Joseph button
end Knife 15cm




Killian Nassal Speculum


35mm
HMI-PS-2001




55mm
HMI-PS-2002



65mm
HMI-PS-2003




90mm
HMI-PS-2004




Cottle Nasal Speculum


20mm
HMI-PS-2005




30mm
HMI-PS-2006




40mm
HMI-PS-2007



50mm
HMI-PS-2008




Let Us know if you could not find what you are looking for



Joseph Nasal
Saw 17cm Angled Left

HMI-PS-2009



Joseph Nasal
Saw 17cm Angled Right

HMI-PS-2010



HMI-PS-2011

Joseph Nasal Saw
Straight 19cm



HMI-PS-2012



Joseph Nasal Saw
Bayonet Left 19cm



HMI-PS-2013

Joseph Nasal Saw
Bayonet Right 19cm

Let Us know if you could not find what you are looking for



Let Us know if you could not find what you are looking for



HMI-PS-2201

Aufricht
Fiberoptic Retractor
45mm Blade 16cm



HMI-PS-2202

Aufricht Retractor
45mm Blade 16cm



HMI-PS-2203

Freeman Flap Retractor
Fiber Optic, 8" 13mm Wide




HMI-PS-2204

Deaver Retractor Fiber Optic



Marx breast retractor,
27cm: 30 x 125mm

HMI-PS-2301



Emory Endoscopic Retractor 1x7" Blade for 10mm
Scope

HMI-PS-2302



Double End Breast Retractor
24cm 38MM X 70MM
HMI-PS-2303



Abdominal Retractor Saldanha
for lipoabdominoplasty 37cm

HMI-PS-2304



Ferreira Fiber Optic Retractor

Length	Width	Code
90mm	24mm	HMI-PS-2401
90mm	30mm	HMI-PS-2402
150mm	36mm	HMI-PS-2403
190mm	36mm	HMI-PS-2404



Tebbets Fiber Optic Retractor

Length	Width	Code
90mm	24mm	HMI-PS-2405
90mm	30mm	HMI-PS-2406
150mm	36mm	HMI-PS-2407
190mm	36mm	HMI-PS-2408



Fiber Optic Cable 2.25 meter with Storz Connector

HMI-PS-2409



C-Retractor C area 50mm Ext 85mm

HMI-PS-2410



**BIGGS MAMMAPlasty
Retractor**

HMI-PS-2501

Without Fiber Optic
30mm




HMI-PS-2502

With FiberOptic



HMI-PS-2503

Ferreira Retractor
Fiber Optic
108mm length,
22mm Wide Blade



Biggs Face Lift Retractors 23cm

HMI-PS-2504

Fiber optic

HMI-PS-2505

Without Fiber
optic



RIGHT HMI-PS-2601

**AGRIS-DINGMAN BREAST
DISSECTOR SET**
36cm



LEFT HMI-PS-2602

HMI-PS-2603

Breast dissector
angled 33cm



HMI-PS-2604

Breast dissector
Spatulated 33cm



HMI-PS-2605

SILVERSTEIN Breast Dissector
blunt 35 cm



HMI-PS-2606

MACCOLLUM DINGMAN
41cm



HMI-PS-2701

Arcus Marginalis
Dissector 24cm
S-Shaped 8mm wide



HMI-PS-2702

Flap Dissector 22.5cm
Curved 6mm wide



HMI-PS-2703

Frontoglabellar Dissector
22.5cm S-Shaped 7mm wide



HMI-PS-2704

Frontotemporal
Dissector 22.5cm
Curved 10mm wide



HMI-PS-2705

Frontotemporal
Dissector 22.5cm
Straight 10mm wide



HMI-PS-2706

Nerve Hook - Curved
Left - 24cm



HMI-PS-2707

Nerve Hook - Curved Right -
24cm



Let Us know if you could not find what you are looking for

HMI-PS-2708

Nerve Retractor -
Curved Left 17cm



HMI-PS-2709

Nerve Retractor -
Curved Right 17cm



HMI-PS-2710

Periosteal Dissector
25cm Curved 9mm wide



HMI-PS-2711

Periosteal Dissector/Elevator
S-Shaped 9mm wide



HMI-PS-2712

Scalp Elevator 25cm
Curved 12mm wide



HMI-PS-2713

Temporal Line Dissector
23.5cm Straight 10mm wide



HMI-PS-2714

Zygomatic Arch Elevator
22.5cm Straight 5mm wide





Collin Mallet 20cm
 30mm - 220 g
 HMI-PS-2801



Hajek Mallet 21.5cm
 27mm - 140 g
 HMI-PS-2802



Lucae Mallet 21.5cm
 19mm - 230 g
 HMI-PS-2803



Mead Mallet 18cm
 26mm - 320 g
 HMI-PS-2804



Plastic Facing Mallet 18cm
 26mm - 320 g
 HMI-PS-2805



Fiber Mallet 25cm
 35mm - 150 g
 HMI-PS-2806

Let Us know if you could not find what you are looking for



HMI-PS-2901

**Asch Septum
Straightening
Forceps 22cm**



Straight
HMI-PS-2902

Cottle Walsham 23cm



Left
HMI-PS-2903

Right
HMI-PS-2904



HMI-PS-2905

**Jansen Middleton
Septum Forceps
Double Action 19cm**



HMI-PS-2906

**TC RUBIN
morselizer
18.5cm**



Cartilage Crusher, 65x30mm

HMI-PS-3001



Cartilage Crusher

HMI-PS-3002




Cartilage Crusher Serrated

HMI-PS-3003



Sheen Cartilage
Crusher Morselizer

HMI-PS-3004



Sheen Grid
10cm x 7cm

HMI-PS-3005



Teflon Cartilage Cutting
Block 85x55x10mm

HMI-PS-3006

HMI-PS-3101

Bipolar Forceps
12 cm . 1mm tip



HMI-PS-3102

Bipolar Forceps
15 cm . 1.5mm tip



HMI-PS-3103

TC Crow Coagulating
Dressing Forceps 23cm



HMI-PS-3104

Bipolar Cable
3m cable length



HMI-PS-3105

Desmaress Lid Retractor
12 insulated



HMI-PS-3106

US Army Retractor
Insulated



HMI-PS-3107

Jaeger Lid Plate Insulated 11cm



Let Us know if you could not find what you are looking for

Thank you for going through Catalog



www.healermedinstruments.com
Email: Healermedins@gmail.com
Whatsapp & Call: +92 3034515151